

{pgslideshow id=3|width=250|}

- The Glove Caddy can hold glove boxes up to 5 1/2" W x 3 1/2" H x 10" L.
- Provides convenient storage for glove box.
- Easy access will increase use of gloves.
- Inexpensive, functional, and easy to install with no drilling or modification to portable machine.
- Original design allows for ease of box insertion, but prevents glove box from dislodging while machine is in motion.
- New Glove Caddy **powder-coat option!** Color matches your GE portable. Looks great on the new AMX700 digital portable!

For all AMX portable x-ray machines

Part# GBH1 16g aluminum \$69.00 [Add to Shopping Cart](#))

Part# GBH1P Powder coated \$69.00 [Add to Shopping Cart](#))

{gallery}glovecaddy{/gallery}